

CT Syringes



A01103/A01114
 Compatible with Medrad VCT610、
 ECT710, For CT
 Content : 1-200ml syringe
 1-Connecting tube
 1-QFT or 1-Spike



A03103/A03114
 Compatible with Medrad Stellant CT
 Content : 1-200ml syringe
 1-Connecting tube
 1-QFT or 1-Spike



A03106/A03107
 Compatible with Medrad Stellant
 Dual CT
 Content : 2-200ml syringes
 1-Y Connecting tube
 2-QFT or 2-Spikes

MAIWEI

Shenzhen Maiwei Biotech Co., Ltd.

No. 2-10, Jinlong Blvd South, Pingshan District, Shenzhen

Tel: +86 755 8522 2663 Fax: +86 755 8522 2785

E-mail: export@maiweibio.com Website : www.maiweibio.com




CT Syringes



B01103/B01114
Compatible with LF CT9000&CT9000ADV
Content : 1-200ml syringe
1-Connecting tube
1-QFT or 1-Spike



B04106/B04107
Compatible with LF OptiVantage Dual CT
Content : 2-200ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes



C01103/C01114
Compatible with Nemoto A-25,A-60,
Dual Shot, For CT
Content : 1-100ml syringe
1-Connecting tube
1-QFT or 1-Spikes

CT Syringes



C02103/C02114
Compatible with Nemoto A-60,
Dual Shot ,For CT
Content : 1-200ml syringe
1-Connecting tube
1-QFT or 1-Spike



C05104/C05106/C05107
Compatible with Nemoto Dual Shot
CT
Content : 1-100ml syringe
1-200ml syringe
1-Y Connecting tube
1-QFT 1-Spike
Or 2-QFT/2-Spikes



C06104/C06106/C06107
Compatible with Nemoto Dual Shot
CT
Content : 2-200ml syringes
1-Y Connecting tube
1-QFT 1-Spike
Or 2-QFT/2-Spikes



CT Syringes



C07106/C07107
Compatible with Nemoto Dual Shot CT
Content : 2-100ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes

CT Syringes



E03106/E03107
Compatible with Imaxeon Salient,
Dual,CT
Content : 2-190ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes



E03103/E03114
Compatible with Imaxeon Salient, CT
Content : 1-190ml syringe
1-Connecting tube
1-QFT or 1-Spike



F01103/F01114
Compatible with EZEM Empower , CT
Content : 1-200ml syringe
1- Connecting tube
1-QFT or 1-Spike



CT Syringes



F02106/F02107
Compatible with EZEM Empower,
CT
Content : 2-200ml syringe
1-Y Connecting tube
2-QFT or 2-Spike



H02103/H02114
Compatible with Seacrown, For CT
Content : 1-200ml syringe
1-Connecting tube
1-QFT or 1-Spike



H02106/H02107
Compatible with Seacrown Dual, For
CT
Content : 2-200ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes

CT Syringes



J01103/J01114
Compatible with Antmed , For CT
Content : 1-200ml syringe
1-Connecting tube
1-QFT or 1-Spike



J01106/J01107
Compatible with Antmed Dual , For
CT
Content : 2-200ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes



N01103/N01114
Compatible with Sino , For CT
Content : 1-200ml syringe
1-Connecting tube
1-QFT or 1-Spike

CT Syringes



N01106/N01107
Compatible with Sino Dual , For CT
Content : 2-200ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes



~~P01103/P01114~~
Compatible with ~~Apostar~~ , For CT
Content : 1-100ml syringe
1-Connecting tube
1-QFT or 1-Spikes



~~P05106/P05107~~
Compatible with ~~Apostar Dual~~ , For CT
Content : 2-100ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes

CT Syringes



P02103/P02114
Compatible with Nemoto A-60,
Dual Shot ,For CT
Content : 1-200ml syringe
1-Connecting tube
1-QFT or 1-Spike



~~P03106/P03107~~
Compatible with ~~Nemoto Dual Shot~~
~~CT~~
Content : 1-100ml syringe
1-200ml syringe
1-Y Connecting tube
2-QFT or 2-Spikes



~~P04106/P04107~~
Compatible with ~~Nemoto Dual Shot~~
CT
Content : 2-200ml syringes
1-Y Connecting tube
2-QFT or 2-Spikes



Angiography Syringes



A02202/A02203
Compatible with Medrad Mark V, DSA
Content : 1-150ml
1-QFT or 1-Spike



A05202/A05203
Compatible with Medrad Mark VII/AT, DSA
Content : 1-150ml syringe
1-QFT or 1-Spike



B02202/B02203
Compatible with LF Angiomat 6000, DSA
Content : 1-150ml syringe
1-QFT or 1-Spike

Angiography Syringes



B03202/B03203
Compatible with LF Angiomat ~~Humera~~
~~DSA~~
Content : 1-150ml syringe
1-QFT or 1-Spike



C02202/C02203
Compatible with Nemoto 120S, DSA
Content : 1-125ml syringe
1-QFT or 1-Spike



C03202/C03203
Compatible with Nemoto Injector, DSA
Content : 1-150ml syringe
1-QFT or 1-Spike



Angiography Syringes



H01202/H01203
Compatible with Seacrown, For DSA
Content : 1-150ml
1-QFT or 1-Spike



N01202/N01203
Compatible with Sino, For DSA
Content : 1-150ml
1-QFT or 1-Spike



P01202/P01203
Compatible with Apostar , For DSA
Content : 1-150ml syringe
1-QFT or 1-Spike



MRI Syringes



A01305
Compatible with Medrad MRI Spectris solaris
Content : 2-65ml syringe
1-Y Connecting tube
2-Spike

MRI Syringes



B01305
Compatible with LF OPTISTAR, MRI
Content : 2-60ml syringe
1-Y Connecting tube
2-Spike



A02305
Compatible with Medrad MRI Spectris solaris
Content : 1-65ml ;
1-115ml syringe
1-Y Connecting tube
2-Spike



C01305
Compatible with Nemoto Sonic Shot, MRI
Content : 2-60ml syringe
1-Y Connecting tube
2-Spike

MRI Syringes



H01305
Compatible with Seacrown, For MRI
Content : 2-65ml syringe
1-Y Connecting tube
2-Spike



F01305
Compatible with EZEM Empower,
MRI
Content : 2-100ml syringe
1-Y or Connecting tube
2-Spike



N01305
Compatible with Sino, For MRI
Content : 2-65ml syringe
1-Y Connecting tube
2-Spike




Pressure Connecting Tube

Pressure Connecting Tube



SOL150
Coil tube,350psi ,150cm with 1 Checkvalve



S0H075/S0H100/S0H120/S0H150/S0H180/
S0H200
Straight tube,1200psi ,braided inside
75cm/100cm/120cm/150cm/180cm/200cm



SOL030
Straight tube,350psi ,30cm with 1 Checkvalve



D0L180
Y Coil tube,350psi, with1
or 2 Check valves, 180cm



D0L150
Y Coil tube,350psi with1 or 2 Check valves
,150cm



Accessories



SOL002
3-way Stopcock
(Male Cap closed and Female Cap open)



SOL001
J-tube/Quick Fill Tube



SOL004
High Pressure 3-way Stopcock
(Male Cap closed and Female Cap open)

